

Sim Pack 1:



Fluorescence

Version 1.0

© Simulators for Teaching, 2012. All rights reserved.



Copyright 2012

www.sim4t.com

Learning Outcomes:

After reading these notes you will learn about the following terms

- **Luminescence (fluorescence & phosphorescence)**
- **Steady state spectroscopy**
- **Fluorescence quantum yield**
- **Fluorescence excited state lifetime**

Steady State Spectrometer

Luminescence spectra (intensity of luminescence as a function of wavelength) can be produced on a steady state spectrometer as shown in Figure 2:




Figure 2: Steady state spectrometer

To produce an emission spectrum the λ_{ex} monochromator (see Figure 2) is held constant and the λ_{em} [or λ_a ($a = \text{analysis}$)] monochromator is scanned.




Figure 3: Steady state emission spectrum

An emission spectrum is therefore the intensity of luminescence (I_L) as function of wavelength (see Figure 3).

Fluorescence Lifetime (Lifetime of the Excited State)

19.

If we use light of intensity I_0 this creates an excited state population $[^1M^*]$ (M absorbs the light to create $^1M^*$), [See Figure 6]

$[^1M^*]$ will remain constant (i.e., in a steady state) until the light is removed

At the instant light is removed the concentration of excited states is $[^1M^*]_0$




Figure 6: decay of excited state population



This document is protected under copyright law. Unauthorised reproduction is prohibited. This applies to all forms of the document including, but not limited to, printed and electronic versions.

© Simulators for Teaching, 2012. All rights reserved



Copyright 2012

www.sim4t.com